Xeo[®]SA

NON-ABLATIVE SKIN REVITALISATION

A customisable, multi-application light-based platform designed to treat the widest range of today's most common non-ablative skin concerns.













xeo SA's comprehensive offering of intense pulsed light options provides a complete range of treatments associated with the signs of photo-aging and other skin imperfections.



LIMELIGHT

Customisable approach for 3-in-1 pigment & vascular treatments

- + Programmable wavelengths from 520 1100 nm: 3 distinct pre-set programs provide excellent hemoglobin and melanin absorption for superior outcomes, uniform treatments, and ease of use.
- + Fast treatments: 10 mm x 30 mm tip provides fast & effective treatments up to 2 Hz.
- + Real time automatic calibration: Seamless monitoring and adjustments ensure accurate output, pulse after pulse.
- + Sun mode: Reduces cooling to the skin, allowing a more aggressive treatment of superficial lentigines while maintaining a high safety margin.







Photos Courtesy of Dr. Jill Lezaic



Photo-rejuvenation market is expected reach \$2.5B by 2025¹



Aging skin and increased skin disorders are driving demand for photorejuvenation¹



Hair removal is in the top 5 non-surgical procedures²



ACUTIP

Targeted and precise energy delivery to discrete lesions

- **+ Optimised wavelength spectrum:** Ranging from 500 635 nm, AcuTip's wavelength range is specifically optimized for hemoglobin and melanin.
- → Increased visibility: AcuTip's small contact tip (6.35 mm) offers exceptional visibility and unobstructed views of the treatment area to trace discrete lesions rapidly and effectively.
- + Ergonomic handpiece: Designed to minimize user fatigue, the AcuTip handpiece is light, ergonomic, and easy to use.



PROWAVE LX

Power of an Alexandrite Laser, Delivered with the Ease of an IPL

- **Predictable outcomes:** ProWave LX provides the clinical efficacy on finer, lighter hair by harnessing high peak powers, short pulse widths and consistent energy monitoring throughout the treatment.
- + Enhanced patient comfort: Aggressive, adjustable contact cooling from 5 - 20 °C for superior comfort, safety, and efficacy on skin types I - IV.
- + Unique handpiece design: Low profile design enhances operator visibility during the treatment, while maintaining fast treatments with the 10×30 mm tip.
- + Intuitive user interface: Designed with 3 preset options, ProWave LX features an intuitive user interface, ideal for delegating and providing best-in-class hair removal on all skin types.



¹Global Market Insights, Skin Rejuvenation Market Size, COVID-19 Impact Analysis, Regional Outlook, Application Potential, Price Trends, Competitive Market Share & Forecast, 2022 - 2028.

²Aesthetic Plastic Surgery National Databank Statistics 2020.

© 2022 Cutera ANZ. All rights reserved.